
CHEMICAL COMPATIBILITY CHART



Chlorosulphonic acid Hydrochloric acid (aqueous) Hydrofluoric acid (aqueous) Hydrogen chloride (anhydrous) Hydrogen fluoride (anhydrous) Nitric acid Oleum Phosphoric acid Sulfuric acid

Acetic acid Butyric acid (n-) Formic acid Propionic acid Rosin Oil Tall oil

Group 3: Caustics Caustic potash solution Caustic soda solution

Aminoethylethanolamine Aniline Diethanolamine Diethylenetriamine Diisopropanolamine Dimethylamine Ethylenediamine Hexamethylenediamine 2-Methyl-5-ethylpyridine Monoethanolamine Monoisopropanolamine Morpholine Pyridine Triethanolamine Triethlamine Triethylenetetramine Trimethylamine

Group 5: Halogenated Compounds

Allyl chloride Carbon tetrachloride Chlorobenzene Chloroform Chlorohydrines, crude Dichlorobenzene (o-) Dichlorobenzene (p-) Dichlorodifluoromethane Dichloroethyl ether Dichloropropane Dichloropropene Ethyl chloride Ethylene dibromide Ethylene dichloride Methyl bromide Methyl chloride Methylene chloride Monochlorodifluoromethane Perchloroethylene Propylene dichloride 1,2,4-Trichlorobenzene 1.1.1-Trichloroethane

Trichloroethylene Trichlorofluoromethane Allyl alcohol Amyl alcohol 1.4-Butanediol Butyl alcohol (iso, n, sec, tert) Butylene alvcol Corn syrup Cyclohexyl alcohol Decyl alcohol (n, iso) Dextrose solution Diacetone alcohol Diethylene glycol Diethylene glycol dimethyl ether Diethylene glycol monobutyl ether Diethylene glycol monoethyl ether Diethylene glycol monomethyl ether Diisobutyl carbitol Dipropylene glycol Dodecanol Ethoxylated dodecanol Ethoxylated pentadecanol Ethoxylated tetradecanol Ethoxylated tridecanol Ethoxytriglycol Ethyl alcohol Ethyl butanol 2-Ethylbutyl alcohol 2-Ethylhexyl alcohol Ethylene glycol Ethyleneglycol monobutyl ether Ethylene glycol monoethyl ether Ethylene glycol monomethyl ether Furfuryl alcohol Glycerine Heptanol Hexanol Hexylene glycol Isoamyl alcohol Isooctyl alcohol Methyoxytriglycol Methyl alcohol Methylamyl alcohol Molasses, all Nonanol Octanol Pentadecanol Polypropylene glycol methyl ether Propyl alcohols (n, iso) Propylene glycol Sorbitol Tetradecanol Tetraethylene glycol Tridecyl alcoho Triethylene glycol Undecanol

Group 6: Alcohols, Glycols and Glycol Ethers Group 7: Aldehydes Acetaldehyde Acrolein (inhibited) Butyraldehyde (n, iso) Crotonaldehyde Decaldehyde (n, iso) 2-Ethyl-3-propylacrolein Formaldehyde solutions Furfural Hexamethylenetetramine Isooctyl aldehyde Methyl butyraldehyde Methyl formal Paraformaldehyde Valeraldehyde Group 8: Ketones Acetone

Acetophenone Camphor oil Cyclohexanone Diisobutyl ketone Isophorone Mesityl oxide Methyl ethyl ketone Methyl isobutyl ketone

Group 9: Saturated Hydrocarbons

Butane Cyclohexane Ethane Heptane Hexane Iso-butane Liquified natural gas Liquified petroleum gas Methane Nonane n- Paraffins Pentane Petrolatum Petroleum ethers Petroleum naphtha Polybutene Propane Propylene butylene polymer

Group 10: Aromatic Hydrocarbons

Benzene Cumene p-Cymene Coal tar oil **Diethylbenzene** Dodecyl benzene Dowtherm Ethylbenzene Naphtha, coal tar Naphthalene (includes molten) Tetrahdyronaphthalene Toluene Triethyl benzene Xylene (m-, o-, p-)

Butylene 1-Decene Dicyclopentadiene Diisobutylene Dipentene Dodecene 1-Dodecene Ethylene Liquified petroleum gas 1-Heptene 1-Hexane Isobutylene Nonene 1-Octene 1-Pentene Polybutene Propylene Propylene butylene polymer Propylene tetramer (dodecene) 1-Tetradecene 1-Tridecene Turpentine 1-Undecene

Group 13: Esters

Amyl acetate

Amyl tallate Butyl acetates (n, iso, sec)

Butyl benzyl phthalate Castor oil Croton oil

Epoxidized vegetable oils Ethyl acetate

Ethylene glycol monoethyl ether acetate Ethylhexyl tallate

14: Monomers and Polymerizable

Ethyl acrylate (inhibited) 2-Ethylhexyl acrylate (inhibited) Isodecyl acrylate (inhibited)

Methyl methacrylate (inhibited) o-Propiolactone Styrene (inhibited)

Vinylidene chloride (inhibited) Vinyl toluene

Isoprene (inhibited) Methyl acrylate (inhibited)

Vinyl acetate (inhibited) Vinyl chloride (inhibited)

Dibutyl phthalate Diethyl carbonate Dimethyl sulfate Dioctyl adipate Dioctyl phthalate

Ethyl diacetate

Glycol diacetate Methyl acetate Methyl amyl acetate Neatsfoot oil

Propyl acetates (n, iso) Resin oil

Acrylic acid (inhibited)

Butadiene (inhibited) Butyl acrylate (n, iso)

Fish oil

Olive oil

Peanut oil

Sperm oil Tallow Tanner's oil Vegetable oil Wax, carnauba

Acrvlonitrile

Carbolic oil

Nonylphenol Phenol

Group 16: Alkylene Oxides

Ethylene Oxide

Propylene Oxide

Acetone cyanohydrin Ethylene cyanohydrin

Ammonium hydroxide

Diethyl ether (ethyl ether)

Cresols

Group 18: Nitriles

Acetonitrile

Adiponitrile

Group 19: Ammonia

Group 20: Halogens

Bromine Chlorine

1, 4-Dioxane

Isopropyl ether Tetrahydrofuran

Creosote, coal tar

Sova bean oil

Asphalt Gasolines CasingeadAutomotive Aviation Jet Fuels JP-1 (kerosene) JP-3 JP-4 JP-5 (kerosene, heavy) Kerosene Mineral spirits Naphtha (non aromatic) Naphtha Solvent Stoddard solvent • VM&P Oils Absorption oil Clarified oil Crude oil Diesel oil Fuel oil No. 1 (kerosene) No. 1-D No. 2 No. 2-D No. 4 No. 5 No. 6

Lubricating oil Mineral oil Mineral seal oil

- Motor oil
- Penetration oil Range oil
- Road oil
- Spindle oil
- Spray oil Transformer oil
- Turbine oil
- Acetic anhvdride Propionic anhydride

Group 23: Sulphur, molten

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any ners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details



					Polyethylene			
	Polypropylene		Polycarbonate		LD		HD	
Chemical	20°C	60°C	20°C	60°C	20°C	60°C	20°C	60°C
Acetic acid solution <90%	R	R	R	R	R	R	R	R
Acetic acid, glacial	R	U	U	U	А	U	R	L
Acetone	R	R	U	U	E	U	R	R
Ammonia	R	R	U	U	R	R	R	R
Aniline	R	R	U	U	U	U	R	R
Benzene	А	U	U	U	U	U	U	U
Boric acid	R	R	R	R	R	R	R	R
Bromine	U	U	U	U	U	U	U	U
Carbon tetrachloride	U	U	U	U	U	U	А	U
Chlorinated solvents	U	U	U	Ŭ	Ű	U	U	Ŭ
Chloroform	U	U	U	U	Ū	U	Ŭ	Ŭ
Cresols	U	Ŭ	Ŭ	Ŭ	Ŭ	Ŭ	Ŭ	Ŭ
Diethyl ether	U	Ŭ	U	Ŭ	Ŭ	U	Ŭ	Ŭ
Diethylene glycol	R	R	R	Ŭ	R	R	R	Ē
Ethers	U	Ŭ	A	Ŭ	E	U	A	U
Ethyl acetate	Ā	A		Ũ	Ā	U	A	Ŭ
Ethanol	R	R	R	R	R	U	R	R
Ethylene glycol	R	R	R	U	R	R	E	E
Formaldehyde solutions 40%	R	R	R	R	R	E	R	R
Formic acid	R	R	R	R	R	E	R	R
Hydrochloric acid conc.	E	E		U	R	R	SW	SW
Hydrofluoric acid 1-60%	E	E	20%	U	R	R	R	300
Hydrofluoric acid conc.	E	E	20%	U	R	U	R	
	R	R	R	R	R		R	
Hydrogen peroxide 30-90%	R	R	п	n		U	R	A
Hypochlorous acid					A	U		R
Magnesium hydroxide	R	R			R	R	R	R
Maleic acid	R	R			R	R	R	R
Mineral oil	R	R	R	R	A	U	R	E
Nitric acid <25%	R	R	R	U	R	R	R	R
Nitric acid 25-70%	R	E	U	U	A	U	L	U
Nitric acid >70	E	U	U	U	U	U	U	U
Oxalic acid	R	A	R	R	R	R	R	R
Paraffin	R	A	R	R	A	U		
Petroleum ether	A	U	R	R	U	U	U	U
Phosphoric acid <30%	R	R	R	R	R	R	R	R
Phosphoric acid 30-85%	R	R	R	R	R	A	R	R
Phosphoric acid >85%	R	R	R	R	A	U	R	А
Potasium hydroxide	50%	10%	U	U	R	R	R	30%
Salicylic acid	R	R	R	R	R	R	R	R
Sodium hydroxide solution	50%	10%	U	U	R	R	R	30%
Sulfuric acid <60%	R	R	R	R	R	R	R	R
Sulfuric acid 60-70%	R	U	R	R	R	А	R	А
Sulfuric acid >70%	E	U	U	U	А	U	А	L
Silicone oil	R	R	R	R	А	U	E	Е
Stearic acid	R	R	R	R	R	R	R	R
Toluene	U	U	U	U	U	U	U	U
Trichloroethylene	U	U	U	U	U	U	V	U
Water	R	R	R	R	R	R	R	R

CHEMICAL RESISTANCE OF PLASTICS

Key

R Resistant

SW Swelling occurs

L Limited resistance

% Max concentration

A Slowly attacked (not recommended for long term storage)E Environmental stress cracking

V Vapour diffusion U Unsuitable

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.